



Floor Scales

Marsden M-420

Portable, High Accuracy Floor Scale



- Highly robust - will outlive most scales
- Large, clear display and easy to use functions
- Suitable for hospitals and health visitors
- Capacity of 220kg, accurate to 50g



The M-420 is Marsden's flagship floor scale. It's a best-seller around the world, thanks to its hard-wearing design and fine gradations of 50g. An optional carry case makes it the perfect portable weighing solution.

Accessories:



CC-420



Bluetooth

Capacity:	220kg
Graduations:	50g<150kg>100g
Dimensions:	320mm x 310mm x 60mm
Power Supply:	Mains or 6 x AA batteries

Adaptor Specification:	12V 1A 2.5mm +ve
Battery Life:	Over 100 hours of continuous use
Weight of Scale:	3.4kg

Marsden M-430

Robust and Portable Floor Scale



- Highly robust - will outlive most scales
- Large, clear display and easy to use functions
- Suitable for GP surgeries and health centres
- Capacity of 220kg, accurate to 200g



Identical in design to our popular M-420, the Marsden M-430's 200g graduations are ideal for GP surgeries and health centres. The reinforced aluminium structure means this scale is lightweight yet built to provide years of precision weighing.

Accessories:



CC-420



Bluetooth

Capacity:	220kg
Graduations:	200g
Dimensions:	320mm x 310mm x 60mm
Power Supply:	Mains or 6 x AA batteries

Adaptor Specification:	12V 1A 2.5mm tip +ve
Battery Life:	Over 100 hours of continuous use
Weight of Scale:	3.4kg

Marsden M-425

Portable Floor Scale with Separate Display



- Highly robust scale - portable and built to last
- Remote handheld weight display for discreet readings
- Suitable for hospitals, school nurses and health visitors
- Easy to use and accurate to 50g



For weight sensitive environments, use the M-425. Based on our rugged, portable M-420 and M-430 floor scales, it's and built to last. With a clear handheld weight display, it's perfect for health visitors, school nurses and family doctors.

Accessories:



CC-420



Bluetooth

Capacity:	220kg	Adaptor Specification:	12V 1A 2.5mm tip +ve
Graduations:	50g<150kg>100g	Battery Life:	Over 100 hours of continuous use
Dimensions:	320mm x 310mm x 60mm	Weight of Scale:	3.4kg
Power Supply:	Mains or 6 x AA batteries		

Marsden M-545

Entry Level Portable Floor Scale



- Robust and portable with integrated carry handle
- Remote handheld weight display for discreet readings
- Measures BMI; suitable for hospitals and GP surgeries
- Affordable, high accuracy scale



With a 250kg capacity, 100g accuracy and rugged build quality, the M-545 is a dependable entry-level floor scale. The Body Surface Area (BSA) function is perfect when calculating drug dosages. The indicator is handheld for weight-sensitive environments.

Accessories:



CC-420

Capacity:	250kg	Adaptor Specification:	12V 1A 2.5mm tip +ve
Graduations:	100g	Battery Life:	180 hours of continuous use
Dimensions:	310mm x 310mm x 90mm	Weight of Scale:	6kg
Power Supply:	Mains x 6 x AA batteries		

Marsden M-510

Premium Portable Floor Scale



- Remote indicator for weight sensitive areas
- Three displays and full numeric keypad
- Suitable for hospitals and weight management clinics
- Integrated carry handle



The M-510 is Marsden's premium floor scale for weight sensitive areas. Three displays show weight, height, plus BMI - which is easy to calculate using the full numeric keypad. The rechargeable battery will last for upto 3000 weigh-ins.

Accessories:



CC-420



Wifi



Bluetooth

Capacity:	250kg
Graduations:	100g
Dimensions:	310mm x 310mm x 90mm
Power Supply:	Mains/Rechargeable

Adaptor Specification:	12V 1A 2.5mm tip +ve
Battery Life:	55 hours of continuous use
Weight of Scale:	6kg

Marsden M-530

High Capacity Bariatric Scale for Plus Size Patients



- Wide, low profile platform
- Remote indicator for discreet weighing
- Suitable for hospitals and weight management clinics
- Capacity of either 300kg or 500kg



A sturdy, extra wide base make the M-530 perfect for weighing larger patients. The indicator is separate, perfect for weight sensitive areas, and three displays show weight, height and BMI. The base has a carry handle for transportation.

Accessories:



CC-530



Wifi



Bluetooth

Capacity:	300kg or 500kg
Graduations:	50g<150kg>100g (300kg) / 100g<300kg>200g (500g)
Dimensions:	600mm x 390mm x 50mm
Power Supply:	Mains/Rechargeable

Adaptor Specification:	12V 1A 2.5mm tip +ve
Battery Life:	55 hours of continuous use
Weight of Scale:	8kg