



Wheelchair Scales

Marsden M-610

Portable Wheelchair Weigh Beams



- Lightweight and easy to set up
- Tare and Preset Tare to deduct weight of wheelchair
- Suitable for weighing most wheelchair types
- Accurate to 100g

Our most portable wheelchair scale.

The M-610 is our most portable wheelchair weighing scale. Perfect for nurses and community health visitors, the beams are lightweight with an integrated handle, and an optional carry case is available. For even greater accuracy, deduct the wheelchair weight using Tare/Preset Tare.

Each beam has a handle and fold-in ramps for portability, and carry cases are available. The M-610 is easy to use, accurate to 100g and perfect for community nurses and healthcare visitors, as well as care homes. Use the full numeric keypad to input the patient's height so BMI can be calculated.

Accessories:



Bluetooth



Wifi



TP-2100



CC-610

Capacity:	300kg	Adaptor Specification:	12V 1A 2.5mm tip +ve
Graduations:	100g	Battery Life:	55 hours of continuous use
Dimensions:	1200mm x 150mm x 95mm	Weight of Scale:	8.5kg per beam
Power Supply:	Mains/Rechargeable		

Marsden M-620

High Capacity Folding Wheelchair Scale



- Large movable weighing platform
- Folds in for easy storage and transportation

- Tare and Preset Tare to deduct weight of wheelchair
- Capacity of 300kg - suitable for most patients

Weigh almost any wheelchair, almost anywhere.

The M-620 is a wheelchair scale with a large spacious weighing platform to accommodate almost every size of wheelchair. For easy transportation and storage the ramp folds in and the platform folds up into a compact unit. Tare/Preset Tare mean you can deduct the weight of the wheelchair to get an accurate weight reading for your patient.

Accessories:



Wifi



TP-2100



Bluetooth



Capacity:	300kg
Graduations:	100g
Dimensions:	700mm x 800mm x 80mm
Power Supply:	Mains/Rechargeable

Adaptor Specification:	12V 1A 2.5mm tip +ve
Battery Life:	55 hours of continuous use
Weight of Scale:	50kg

Marsden M-650 Range

Professional Wheelchair Scales



- Customisable wheelchair scale
- One piece platform with gentle incline
- Tare and Preset Tare to deduct wheelchair weight
- Accurate to 100g

Tailor this scale to your needs.

The M-650 has a large, low profile platform that will accommodate almost any wheelchair. Our intuitive indicator has a full numeric keypad for easy BMI calculation. For easier use, the platform is available with an indicator column (M-651) so that the weight display is easy to read at waist height.

For patients able to stand, two handrails (M-652) make weighing a struggle-free process, with the weight display mounted on one handrail. Alternatively, the M-653 has one hinged handrail.

The scale provides accurate weight readings to 100g and has a 300kg capacity. BMI can be calculated by inputting the patient's height. The indicator is powered by an internal rechargeable battery for extra portability.

For patients with visual-spatial issues, choose the all-white model.

Accessories:



Wifi



TP-2100



Bluetooth



M-650



M-651



M-652



M-653

Capacity:	300kg
Graduations:	100g
Overall Dimensions:	1150mm x 800mm x 66mm
Weighing Area:	950mm x 740mm

Power Supply:	Mains/Rechargeable
Adaptor Specification:	12V 1A 2.5mm tip +ve
Battery Life:	55 hours of continuous use
Weight of Scale:	28.6kg