



## Column Scales

# Marsden M-100

Column Scale with Integrated Height Measure



- Unique corner design
- Large platform - easier to step onto
- Integrated height rod - a complete measuring station
- Suitable for hospitals and weight management programmes

## High capacity weighing and measuring.

Weighing and measuring a patient become one process with the M-100: Its integrated height measure and large base - suitable for plus size patients - mean getting the readings you need is now easier. BMI is easy to calculate, too, with a full numeric keypad for entering the patient's height. Better still, the M-100 has been uniquely designed to allow it to fit neatly into the corner of a room, saving space. For greater convenience, an automatic height measure and printer are available for this scale.



### Accessories:



Bluetooth



Wifi



TP-2100



HM-201D

Capacity: 300kg  
 Graduations: 50g < 150kg > 100g  
 Dimensions: 1240mm x 440mm x 460mm  
 Power Supply: Mains/Rechargeable

Adaptor Specification: 12V 1A 2.5mm tip +ve  
 Battery Life: 55 hours of continuous use  
 Weight of Scale: 15kg

# Marsden M-110

Easy To Use Column Scale with BMI



- High accuracy, easy to use column scale
- Full numeric keypad for quick BMI calculation
- Optional height measures and printer
- Capacity of 250kg

## Simple weighing to complete measuring station.

The M-110 is a premium quality column scale ideal for use in clinics and on hospital wards. It's easy to operate and a rechargeable battery means it can be used anywhere.

Wheels at the bottom of the column enable safe and easy transportation between rooms.

For extra convenience choose the optional height measures - the HM-201M manual height measure, or the HM-201D, which automatically sends height data to the scale.

These height measures, along with optional wireless connectivity and printer, make the M-110 a complete measuring station.



### Accessories:



Bluetooth



TP-2100



Wifi



HM-201D



HM-201M

|               |                       |
|---------------|-----------------------|
| Capacity:     | 250kg                 |
| Graduations:  | 100g                  |
| Dimensions:   | 935mm x 310mm x 310mm |
| Power Supply: | Mains/Rechargeable    |

|                        |                            |
|------------------------|----------------------------|
| Adaptor Specification: | 12V 1A 2.5mm tip +ve       |
| Battery Life:          | 55 hours of continuous use |
| Weight of Scale:       | 9kg                        |

# Marsden M-125

Entry Level Digital Column Scale



- Robust scale for hospital wards and GP surgeries
- Large, clear 3cm digital display
- BSA and BMI calculation
- Accurate to 100g

## Entry level price, Class III accuracy.

The M-125 will help you fulfil all basic medical weighing needs - with 100g accuracy, 250kg capacity and simple weighing functions. Use the scale to calculate a patient's BMI or Body Surface Area (BSA) - a feature unique to Marsden scales. Results are displayed on the large, clear LCD display.

As Marsden's most affordable column scale, the M-125 is ideal for GP surgeries, health centres and clinics. The BSA can benefit renal, oncology and other areas where clear BMI, BSA and weight information is needed. Outpatient departments can benefit from the large display and sturdy construction.



### Accessories:



HM-201M



|               |                           |
|---------------|---------------------------|
| Capacity:     | 250kg                     |
| Graduations:  | 100g                      |
| Dimensions:   | 960mm x 310mm x 310mm     |
| Power Supply: | Mains or 6 x AA batteries |

|                        |                             |
|------------------------|-----------------------------|
| Adaptor Specification: | 12V 1A 2.5mm tip +ve        |
| Battery Life:          | 180 hours of continuous use |
| Weight of Scale:       | 9kg                         |