

AQUATHERM

Heating mattress

NIDO

Baby comfort program



Aquatherm

A soft and warm nest for the little patient to make him feel like in the maternal womb...

When the baby is born, also the premature one, it requires different treatments destined to facilitate the growth phase and especially to avoid trauma that could affect the physical state of the baby.

At this stage it is very helpful to use "Aquatherm" because it provides an additional beneficial source of heat to the little patient.

In addition, the Baby Bundle, flexible and easily shapable, is a gentle and effective aid to the newborn to find his natural position.

The "Aquatherm" system is portable, lightweight and easily positionable:

The Control Unit is equipped with feet that make it stable on any flat surface as well as a support to place it on vertical surfaces.

The Mattress can be placed on a bench care (baby changing table) as well as inside an incubator or wherever you need a heated flat surface. The "Aquatherm" system, by virtue of its property to store heat, can be used with the baby crib "Nido" during the transport between different wards within the hospital.



Control Unit

Unit that proportionally servo-controls the rate of heat supplied to the mattress in order to reach and maintain the temperature set by the operator. It is provided with controls and alarms to ensure at the same time easiness of use and maximum reliability of the system.



Heating mattress

The mattress is the result of a new technology "Aqualgel", that offers all the advantages of a soft and impressed surface around the newborn. It is effective in minimizing the vibration and heat loss towards the surface and able to store the heat even when the control unit is not powered or in case of internal transport between wards.



Nido



The baby crib “Nido” has been studied especially for the comfort and the safety of the newborn and to offer the most practical and functional solution to the health operator. The trolley carrying the baby crib has a solid and safe structure and is made of epoxy painted steel. The baby crib “Nido” has been designed to allow the mother to stay near his child allowing her to remain comfortably relaxed in bed. The rounded shape of the hood aims to guarantee the safety of the newborn and of the nurse while its non-toxic material grants an easy cleaning.

The exclusive cover system guarantees an advanced comfort and the best protection also during transports within the hospital. The standard “Nido” includes a trolley equipped with four castors (two of which are with brakes), a mattress made of fire retardant material and a plastic bassinet that can be disinfected. More than a baby crib, “Nido” has been created to be a real neonatal care system, offering a wide range of accessories.

The bassinet can be placed in Trendelenburg and Fowler position in the simplest way while the handles of the baby crib hood blocking system, positioned at the baby crib’s sides, allow a free and complete rotation of the hood permitting a complete access to the newborn. The structure in transparent plastic material also allows the newborn to receive phototherapy treatment without the need of being moved. The complete freedom of movement is granted by the innovative handle in the frontal position that allows the observation of the patient during transports.





Nido

Accessories



Baby crib hood (P/N 11603A70)



Baby Bundle: a soft and flexible aid, which can be mounted around the baby to help him finding his natural position
(P/N 10828A70) (100x8cm)
(P/N 10829A70) (80x5cm)



Drawers (P/N 11899A70)



Basket (P/N 11895A70)



ID Label Holder "Baby Flag" 12 pcs/set (P/N 11987A70)



Screen for Nido (P/N 12455A70)

Nido & Aquatherm

Aquatherm

Technical Specifications

Power Supply	230 Vac, 50–60Hz
Maximum Power	140 VA
Safety Class	I B
Total weight	4.2 Kg
Dimensions	
Control Unit	19 x 15 x 19 cm
Heating Mattress	35 x 61 x 3 cm

Nido

Technical Specifications

Tilting (Trendelemburg and Fowler)	+/- 12°
Max. Load	10 Kg
Dimensions	80 x 50 x 92 cm
WEIGHT	
trolley + cradle	12 kg
hood	3 kg
Castors:	4
Diameter:	8 cm
Number of breaks:	2



GINEVRI

neo-tech for life

GINEVRI srl
Via Cancelliera, 25/b
00041 Albano Laziale (Rome) – Italy
Tel.: +39 06 93459 330
e-mail: export@ginevri.com
www.ginevri.com

Quality System

ISO 9001:20015
ISO 13485:2016



Authorized Distributor



The specifications in this catalogue are indicative. The company GINEVRI srl reserves the right to make changes, without further notice, to the products described within this catalogue in order to improve reliability, function or design.