



Chair Scales

NEW Marsden M-250

Chair Scale with Lift Assist



- Simple process for patients who need assistance
- Smooth-lifting seat - less effort required by staff and reduced risk of back injuries
- Ergonomically designed chair seat for patient comfort
- Enhances manual handling best practices

Reduced strain on the nurse or carer.

The powered seat of the brand new M-250 could reduce effort needed to lift a patient from sitting to standing by as much as 70%, according to the results of our development and testing process.

The M-250 has been designed to reduce stress and pressure on patients and staff. It makes weighing patients who are unable to stand unaided easier: the sculpted seat lowers the patient to sitting position, and then raises them to standing position.

This, in turn, could reduce staff time - and potentially even the need for two staff members - and the risk of back injuries, caused by patient lifting.

The lifting mechanism is easy to operate via the handheld control and features such as hinged armrests and footrests and precision accuracy mean a better experience for the patient and less stress for staff.



Scan me for more information

Capacity:	250kg
Graduations:	100g
Overall Dimensions:	1085mm x 670mm x 1015mm
Seat Dimensions	440mm x 530mm

Power Supply:	Mains/Rechargeable
Adaptor Specification:	12V 1A 2.5mm tip +ve
Battery Life:	55 hours of continuous use

Marsden M-200

High Capacity, Plus Size Chair Scale



- Comfortable, extra wide seat for larger patients
- Three clear displays for weight, height, BMI

- Hinged armrests and footrests for easy access
- Accurate to 50g



The M-200 is perfect for plus size patients with its 250kg or 300kg capacity options and extra wide, one-piece seat. Hinged armrests and footrests mean easy access and greater patient comfort and safety.

Accessories:



Wifi



TP-2100



Bluetooth

Capacity:	250kg or 300kg
Graduations:	50g<150kg>100g
Seat Dimensions:	530mm x 460mm
Power Supply:	Mains/Rechargeable

Adaptor Specification:	12V 1A 2.5mm tip +ve
Battery Life:	55 hours of continuous use
Weight of Scale:	27kg

Marsden M-210

Professional, Easy To Use Chair Scale



- Easy to use scale for patients with limited mobility
- Three clear displays for weight, height, BMI

- Hinged armrests and footrests for easy access
- Capacity of 250kg



A 250kg capacity, 50g accuracy and an easy to use, split-display indicator are the key benefits of the Marsden M-210. Add to that hinged armrests and footrests to make it easier for patients to get into the seat and you have a versatile chair scale for hospitals and care homes.

Accessories:



Wifi



TP-2100



Bluetooth

Capacity:	250kg
Graduations:	50g<150kg>100g
Seat Dimensions:	400mm x 380mm
Power Supply:	Mains/Rechargeable

Adaptor Specification:	12V 1A 2.5mm tip +ve
Battery Life:	55 hours of continuous use
Weight of Scale:	18kg

Marsden M-225

Entry Level Chair Scale



- Entry level solution for weighing whilst seated
- Large, clear 3cm digital display
- Hinged armrests and footrests for easy access
- Accurate to 100g



The M-225 has all the essential functions you need when weighing patients with limited mobility. A 250kg capacity, 100g accuracy and long battery life make it a convenient solution for weighing most patients. Hinged armrests and footrests ensure easy access.



Capacity:	250kg	Adaptor Specification:	12V 1A 2.5mm tip +ve
Graduations:	100g	Battery Life:	180 hours of continuous use
Seat Dimensions:	400mm x 380mm	Weight of Scale:	18kg
Power Supply:	Mains or 6 x AA batteries		

Marsden M-230

Entry Level Chair Scale with Wider Seat



- Chair scale with wider than standard seat
- Large, clear 3cm digital display
- Hinged armrests and footrests for easy access
- Easy to use for simple, accurate results



The Marsden M-230 shares many similarities of the M-225, with BSA and BMI calculation and accuracy to the nearest 100g. But with a seat that's 60mm wider, this chair scale can accommodate larger patients.



Capacity:	250kg	Adaptor Specification:	12V 1A 2.5mm tip +ve
Graduations:	100g	Battery Life:	180 hours of continuous use
Seat Dimensions:	460mm x 430mm	Weight of Scale:	18kg
Power Supply:	Mains or 6 x AA batteries		