



Fetal Doppler

DIXION Overtone 6000

Manual echograph Overtone 6000 meets the routine examination requirements of obstetricians, midwives and pregnant women. Overtone 6000 is a high Performance model with digital LCD display which offers FHR monitoring in 3 modes: real-time, averaged or manual.

Features:

- ◆ LCD display which offers FHR monitoring in 3 modes: real-time, averaged or manual
- ◆ Ergonomic design, compact and light
- ◆ Built-in rechargeable battery (250 examinations)
- ◆ Accurate FHR detection on LCD display
- ◆ Automatic and manual FHR counting
- ◆ Battery charge indicator
- ◆ Interchangeable ultrasound probes
- ◆ Built-in dynamic loudspeaker / output for recorder or headset
- ◆ Automatic power-off

Specification:

Display	45 mm x 25 mm
Functional modes	Continuous wave ultrasound
Frequency	2 MHz
Transmitter power	< 10 MW/cm ² (Isata)
Heart rate range	50 – 210 bpm
Audio output power	1 W
Work temperature	5 – 40 °C
Power supply	Battery 9 W, continuous work > 4 h
Size	138 mm x 85 mm x 32 mm

** Please contact our company representatives for detailed CE information.*