



Emergency Suction Pump

DIXION Vacus 7305

The model Vacus 7305 is a trolley-style high vacuum and high flow medical suction unit. Vacus 7305, which absorbs the mucus rapidly from the body-cavities or wounds, is widely used in the in-patient department, operation rooms in particular.

Features:

- ◆ The oil-less vacuum source requires no maintenance
- ◆ The capacious storage bottles made of high quality plastics and equipped with a joint ring are easy to access, use and clean
- ◆ The overflow protection device prevents the liquid from entering the pump
- ◆ The hoses, made of non-toxic glass like polymer, are easy to use and clean
- ◆ The air filter reduces contamination to the environment
- ◆ The device is controlled by the regulator and the footswitch
- ◆ Low noise
- ◆ Safe and reliable
- ◆ Long-term use

Specification:

Adjustable vacuum range	0,02 MPa – 0,09 MPa (150 – 680 mmHg)
Flow rate	40 ± 2 l/min
Noise	< 60 dB
Storage bottle	2500 ml x 2 glas pcs or 2000 x 2 plastic pcs
Power supply	220 – 230 V, 50 – 60 Hz
Input	250 VA
Weight	19 kg
Size	320 x 330 x 700 mm